

---

## KALADRUG-R: Laboratory SOP #8

### SSG susceptibility testing *L. donovani* intracellular amastigotes

April 6<sup>th</sup> 2010

Adapted: October 26<sup>th</sup> 2009

ITMA

---



#### A. Material required

1. x 10 mL syringes (x=no. of mice)
2. x 21 ½ gauge needles (green)
3. x 26 3/8 gauge needles (brown)
4. 50 mL universals RT
5. 50 mL universals on ice
6. Forceps and scissors
7. 7 mL and 15 ml Sterilin tubes for drug solutions
8. filter syringes 22µM, suitable for small volumes
9. cavity slides (Lab Tek Chamber Slide System 16-well Nunc Cat. No. 178599, VWR)
10. Crystalline SSG (Albert David); NOT the one used for treating patients
11. 2% starch solution

#### B. Solutions required

1. RPMI 1640 (+ L-glutamine, + 25 mM HEPES) (PAA, Catalog no. E15-842)(=Basic RPMI)
2. Gentamycin 10 mg/mL (Sigma, Cat. No. G1272)
3. penicillin/streptomycin (100x) PAA laboratories Cat. No. P11-010 (contents: 10.000 u Pen/ 10mg/mL Strep)
4. RPMI 1640 + HIFCS 10%v/v + 1%v/v Pen/Strep + 0.5% Gentamycin v/v (=Full RMPI)
5. Phosphate Buffer Saline (PBS)

#### C. Preparation of SSG solutions

!!! SSG is extremely toxic and should be handled with utmost care, if any spilling, it should be cleaned thoroughly with water. Upon inhalation of the powder, seek immediate medical help.

!!! All SbV -solutions should be stored in light-protected containers

Molecular Formula SSG =  $C_{12}H_{17}O_{17}Sb_2.3Na$

Total molecular weight: 745.725 g/mol including 243.5 g SbV

For 100% pure SSG: 1 g SSG contains  $243.5/745.725 \text{ g} = 0.3265 \text{ g SbV}$

For x% pure SSG: 1 g SSG contains  $(x\% * 0.3265) \text{ g SbV present}$

SbV content is batch dependent and should be checked with supplier.

for current batch ITMA: 1g SSG contains 0.316 g SbV (= 97.6%)

1. weigh 0.189873 g SSG ( $\approx 0.060 \text{ g SbV}$  or 60 mg SbV)
2. dissolve in 10 ml PBS
3. incubate for 15-20 min @ 37°C to completely solubilise the powder, the resulting solution should be completely transparent
4. filter this solution over 0.22  $\mu\text{M}$  syringe filter for small volumes
5. Aliquot this stock solution 6 mg/mL SbV in volumes of 1 mL and store these aliquots @ -20°C for max. 3 months = **STOCK SOLUTION A (SbV 6 mg/mL)**
6. When starting a susceptibility assay, thaw 1 stock solution A & remember to manipulate sterile
7. Add 500  $\mu\text{L}$  stock A to 500  $\mu\text{L}$  sterile PBS = **stock solution B (SbV 3 mg/mL)**
8. Add 300  $\mu\text{L}$  stock B to 600  $\mu\text{L}$  sterile PBS = **stock solution C (SbV 1 mg/mL)**
9. Add 100  $\mu\text{L}$  stock B to 900  $\mu\text{L}$  sterile PBS = **stock solution D (SbV 0.3 mg/mL)**
10. The stock solutions A, B, C, D should be FRESHLY diluted 1/100 in FULL RPMI just before adding to *Leishmania* infected macrophages resulting in final concentrations of
  - solution A = 60  $\mu\text{g/mL}$
  - solution B = 30  $\mu\text{g/mL}$
  - solution C = 10  $\mu\text{g/mL}$
  - solution D = 3  $\mu\text{g/mL}$

#### D. Extraction and plating murine peritoneal macrophages

1. Inject 2% starch solution in peritoneum of mice to stimulate macrophage production 48hrs before harvest
2. Prepare sufficient (dependent on no. mouse) basic RPMI 1640 with 1% Gentamycin (50 mL basic RPMI + 0.5 mL Gentamycin) in Falcon tube
3. Distribute prepared medium in syringes, 10 mL/ syringe for normal mice, 6 mL/ syringe for BALB/C mice
4. Collect the mice, spread the mouse on their backs and spray belly with 70% ethanol.
5. Snip mouse where the belly button would be (V. small nick!), and pull off the skin (first pull apart, and then sideways = deglove)
6. Spray the inside with 70% ethanol
7. Put a 21 ½ gauge needle onto a charged syringe

8. Using forceps, pull up the membrane at the sternum and insert the needle (just under membrane, bevel facing upwards)
9. Gently fill the cavity with the RPMI (membrane should be taught)
10. Hold the mouse by base of tail and shoulders and wiggle for 20 secs (allowing the RPMI to rinse the cavity)
11. Face the mouse towards you on its back and insert a 26 3/8 gauge (brown) needle (with the empty syringe) into the mouse's side (bevel facing body)
12. Pull the skin taught against the needle and remove the fluid from the cavity, change sides when one side has been depleted
13. When all fluid has been extracted (6-8 mL), transfer the fluid into a 50 mL Falcon on ice.
14. Spin down the collected macrophages @ 1700 rpm for 10 min @ 4°C (a reddish pellet should result, red indicates the presence of RBC)
15. Remove the SN and add 10 mL full RPMI
16. Resuspend the pellet with plastic pasteur pipette and take aliquot to count (macrophages are slightly larger than RBC's and more granular looking, without characteristic shape of RBC)
17. Adjust concentration to  $5 \times 10^5$  cells/mL using full RPMI.
18. Put 10 mL 100  $\mu$ L/well for cavity, prepare 3 separate cavity slides, one for check @ 24hrs (4 wells/strain), one for untreated control infection 5 days (4 wells/strain) & 1 for treated infection 5 days(4 wells/conc. drug/strain)
19. Incubate @ 37°C, 5% CO<sub>2</sub> for 24 hrs (cells will adhere to surface, not multiply!)

Note: The macrophages can be kept for up to a week after plating if necessary, but one must keep the media covering them fresh and suitably high. They are best used immediately after extraction.

### **E. *In vitro* SSG susceptibility testing of *L. donovani* in peritoneal macrophages**

Control strain to be included for all susceptibility tests:

For *L. donovani* Nepal: BPK282/0 clone 4 or BPK206/0 clone 10

For *L. donovani* India: AG83

1. Count late stationary promastigotes in culture
2. Dilute to  $5 \times 10^6$  / mL for infection at 1 mf: 10 parasite ratio using full RPMI
3. Add 100  $\mu$ L parasite suspension of  $5 \times 10^6$  / mL in each well of the cavity slides.
4. Put infected cultures at 37°C, 5% CO<sub>2</sub>
5. @ 24hrs:
  - For controls slide 24 hrs:
    - Remove medium with pastette
    - Wash with basic RPMI 1640 up to 3 times (fill with pastette)
    - Remove the wells from the slide
    - Methanol fix for about 1 min.

- Shake off the methanol
  - Remove the plastic rings of the slide
  - Giemsa stain (1:10 dilution) for 10 min  
(use 10 mL syringe, 4.5 mL H<sub>2</sub>O + 0.5 mL Giemsa, and drip out over slide)
  - Rinse off Giemsa stain with water
  - Count 100 macrophages and determine % infection
  - For other slides:
    - Prepare 1/100 dilutions of the 4 stocks of drug-solutions in full RPMI (see point C)
    - Cover with 200 µL of WORKSOLUTION drug, 4 wells for each concentration
    - Put back at 37°C, 5% CO<sub>2</sub>
6. @ 5 days: stop infection and check cavity slides as described in step 5

## F. Analysis results

1. Microscopically evaluate infection levels of all control wells and all treated wells: count % infected macrophages and average no. of amastigotes/macrophage.
2. Determine infection index = no. of amastigotes/mf \* % infected mf
3. Determine % survival ± stdev at each drug concentration by infection index in treated wells with infection index of untreated wells
4. Use sigmoid regression analysis to determine corresponding ED50 ± 95% confidence intervals  
SPSS, XLfit, Graphpad prism are software packages that can do this type of non-linear regression analysis.